COPPER IUD INSERTION KIT

This insertion kit contains:

- 1 Pozzi tenaculum forceps
- 1 Cheron forceps
- 1 string scissors

Insertion instructions

The IUD insertion kit should only be used by a qualified healthcare professional.

After inserting the speculum, cleanse and disinfect the cervix using the Cheron forceps.

If necessary, grip the cervix with the Pozzi tenaculum forceps to straighten the uterus and position it in line with the vagina.

Use the hystrometer to estimate the uterine depth.

Insert the IUD according to the manufacturer's instructions included in the IUD packaging.

Once the IUD is inserted, cut the strings to the desired length using the string scissors.

After use

After the procedure, destroy the products used and the packaging according to the applicable procedure at your facility for medical waste management.

Recommendations :

Read the IUD manufacturer's insertion instructions.

Verify the integrity of the packaging as well as the expiry date before use.

These devices are for single use.

They should not be used if the packaging has been opened or damaged, and they cannot be reused. Failure to follow this rule may impair the performance of these devices as well as increase the risk of contamination.

Storage conditions:

Sterile unless the package is opened or damaged. Instruments should not be reused and should be destroyed after use.

Store in a dry cool place, away from humidity and direct sunlight.

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EC REP Obelis S.A, Boulevard Général Wahis 53, 1030 Brussels, Belgium