



SMB Code	U.S.P. Size	Length Inches	Length Metric	Needle Length	Needle Curvature and Profile
SA 6013	6-0	30"	76 cm	13mm	1/2 Circle Round Body (Taper Point)
SA 6016	6-0	30"	76 cm	16mm	1/2 Circle Round Body (Taper Point)
SA 6020	6-0	30"	76 cm	20mm	1/2 Circle Round Body (Taper Point)
SB 6008	6-0	30"	76 cm	8mm	3/8 Circle Round Body (Taper Point)
SB 6010	6-0	30"	76 cm	10mm	3/8 Circle Round Body (Taper Point)
SB 6012	6-0	30"	76 cm	12mm	3/8 Circle Round Body (Taper Point)
SB 6013	6-0	30"	76 cm	13mm	3/8 Circle Round Body (Taper Point)
SB 6016	6-0	30"	76 cm	16mm	3/8 Circle Round Body (Taper Point)
SC 6013	6-0	30"	76 cm	13mm	3/8 Circle Cutting
SC 6016	6-0	30"	76 cm	16mm	3/8 Circle Cutting
SD 6008	6-0	30"	76 cm	8mm	3/8 Circle Reverse Cutting
SD 6010	6-0	30"	76 cm	10mm	3/8 Circle Reverse Cutting
SD 6012	6-0	30"	76 cm	12mm	3/8 Circle Reverse Cutting
SD 6013	6-0	30"	76 cm	13mm	3/8 Circle Reverse Cutting
SE 6017	6-0	30"	76 cm	17mm	1/2 or 3/8 Circle Taper Cut
SA 5013	5-0	30"	76 cm	13mm	1/2 Circle Round Body (Taper Point)
SA 5016	5-0	30"	76 cm	16mm	1/2 Circle Round Body (Taper Point)
SA 5020	5-0	30"	76 cm	20mm	1/2 Circle Round Body (Taper Point)
SSA 5020	5-0	30"	76 cm	20mm	Straight Round Body
SB 5008	5-0	30"	76 cm	8mm	3/8 Circle Round Body (Taper Point)
SB 5010	5-0	30"	76 cm	10mm	3/8 Circle Round Body (Taper Point)
SB 5012	5-0	30"	76 cm	12mm	3/8 Circle Round Body (Taper Point)
SB 5013	5-0	30"	76 cm	13mm	3/8 Circle Round Body (Taper Point)
SB 5016	5-0	30"	76 cm	16mm	3/8 Circle Round Body (Taper Point)
SB 5016DA	5-0	30"	76 cm	16mm	3/8 Circle Round Body (Taper Point) Double Arm
SB 5020	5-0	30"	76 cm	20mm	3/8 Circle Round Body (Taper Point)
SB 5020DA	5-0	30"	76 cm	20mm	3/8 Circle Round Body (Taper Point) Double Arm
SC 5013	5-0	30"	76 cm	13mm	3/8 Circle Cutting
SC 5013DA	5-0	30"	76 cm	13mm	3/8 Circle Cutting Double Arm
SC 5016	5-0	30"	76 cm	16mm	3/8 Circle Cutting
SC 5016DA	5-0	30"	76 cm	16mm	3/8 Circle Cutting Double Arm
SCSP 5012	5-0	30"	76 cm	12mm	3/8 Circle Spatulated Micro Needle
SD 5008	5-0	30"	76 cm	8mm	3/8 Circle Reverse Cutting
SD 5010	5-0	30"	76 cm	10mm	3/8 Circle Reverse Cutting
SD 5012	5-0	30"	76 cm	12mm	3/8 Circle Reverse Cutting
SD 5013	5-0	30"	76 cm	13mm	3/8 Circle Reverse Cutting
SD 5013 DA	5-0	30"	76 cm	13mm	3/8 Circle Reverse Cutting Double Arm
SD 5016	5-0	30"	76 cm	16mm	3/8 Circle Reverse Cutting
SE 5017	5-0	30"	76 cm	17mm	1/2 or 3/8 Circle Taper Cut

SMB Code	U.S.P. Size	Length Inches	Length Metric	Needle Length	Needle Curvature and Profile
SA 4013	4-0	30"	76 cm	13mm	1/2 Circle Round Body (Taper Point)
SA 4016	4-0	30"	76 cm	16mm	1/2 Circle Round Body (Taper Point)
SA 4020	4-0	30"	76 cm	20mm	1/2 Circle Round Body (Taper Point)
SA 4025	4-0	30"	76 cm	25mm	1/2 Circle Round Body (Taper Point)
SA 4026	4-0	30"	76 cm	26mm	1/2 Circle Round Body (Taper Point)
SB 4010	4-0	30"	76 cm	10mm	3/8 Circle Round Body (Taper Point)
SB 4012	4-0	30"	76 cm	12mm	3/8 Circle Round Body (Taper Point)
SB 4013	4-0	30"	76 cm	13mm	3/8 Circle Round Body (Taper Point)
SB 4016	4-0	30"	76 cm	16mm	3/8 Circle Round Body (Taper Point)
SB 4020	4-0	30"	76 cm	20mm	3/8 Circle Round Body (Taper Point)
SC 4013	4-0	30"	76 cm	13mm	3/8 Circle Cutting
SC 4013DA	4-0	30"	76 cm	13mm	3/8 Circle Cutting Double Arm
SC 4016	4-0	30"	76 cm	16mm	3/8 Circle Cutting
SC 4016DA	4-0	30"	76 cm	16mm	3/8 Circle Cutting Double Arm
SC 4022	4-0	30"	76 cm	22mm	3/8 Circle Cutting
SC 4026	4-0	30"	76 cm	26mm	3/8 Circle Cutting
SD 4012	4-0	30"	76 cm	12mm	3/8 Circle Reverse Cutting
SD 4013	4-0	30"	76 cm	13mm	3/8 Circle Reverse Cutting
SD 4013DA	4-0	30"	76 cm	13mm	3/8 Circle Reverse Cutting Double Arm
SD 4016	4-0	30"	76 cm	16mm	3/8 Circle Reverse Cutting
SD 4016DA	4-0	30"	76 cm	16mm	3/8 Circle Reverse Cutting Double Arm
SD 4026	4-0	30"	76 cm	26mm	3/8 Circle Reverse Cutting
SE 4017	4-0	30"	76 cm	17mm	1/2 or 3/8 Circle Taper Cut
SA 3013	3-0	30"	76cm	13mm	1/2 Circle Round Body (Taper Point)
SA 3016	3-0	30"	76cm	16mm	1/2 Circle Round Body (Taper Point)
SA 3020	3-0	30"	76cm	20mm	1/2 Circle Round Body (Taper Point)
SA 3025	3-0	30"	76cm	25mm	1/2 Circle Round Body (Taper Point)
SA 3026	3-0	30"	76cm	26mm	1/2 Circle Round Body (Taper Point)
SA 3030	3-0	30"	76cm	30mm	1/2 Circle Round Body (Taper Point)
SA 3035	3-0	30"	76cm	35mm	1/2 Circle Round Body (Taper Point)
SB 3012	3-0	30"	76cm	12mm	3/8 Circle Round Body (Taper Point)
SB 3013	3-0	30"	76cm	13mm	3/8 Circle Round Body (Taper Point)
SB 3016	3-0	30"	76cm	16mm	3/8 Circle Round Body (Taper Point)
SB 3020	3-0	30"	76cm	20mm	3/8 Circle Round Body (Taper Point)
SB 3030	3-0	30"	76cm	30mm	3/8 Circle Round Body (Taper Point)
SB 3045	3-0	30"	76cm	45mm	3/8 Circle Round Body (Taper Point)
SC 3013	3-0	30"	76cm	13mm	3/8 Circle Cutting
SC 3013DA	3-0	30"	76cm	13mm	3/8 Circle Cutting Double Arm
SC 3016	3-0	30"	76cm	16mm	3/8 Circle Cutting



SMB Code	U.S.P. Size	Length Inches	Length Metric	Needle Length	Needle Curvature and Profile
SC 3016DA	3-0	30"	76cm	16mm	3/8 Circle Cutting Double Arm
SC 3022	3-0	30"	76cm	22mm	3/8 Circle Cutting
SC 3026	3-0	30"	76cm	26mm	3/8 Circle Cutting
SC 3030	3-0	30"	76cm	30mm	3/8 Circle Cutting
SC 3045	3-0	30"	76cm	45mm	3/8 Circle Cutting
SD 3016	3-0	30"	76cm	16mm	3/8 Circle Reverse Cutting
SD 3026	3-0	30"	76cm	26mm	3/8 Circle Reverse Cutting
SE 3017	3-0	30"	76cm	17mm	1/2 or 3/8 Circle Taper Cut
SE 3025	3-0	30"	76cm	25mm	1/2 or 3/8 Circle Taper Cut
SSC 3040	3-0	30"/ 36"	76/90cm	40mm	Straight cutting
SSC 3050	3-0	30"/ 36"	76/90cm	50mm	Straight cutting
SA 2013	2-0	30"	76 cm	13mm	1/2 Circle Round Body (Taper Point)
SA 2016	2-0	30"	76 cm	16mm	1/2 Circle Round Body (Taper Point)
SA 2020	2-0	30"	76 cm	20mm	1/2 Circle Round Body (Taper Point)
SA 2025	2-0	30"	76 cm	25mm	1/2 Circle Round Body (Taper Point)
SA 2026	2-0	30"	76 cm	26mm	1/2 Circle Round Body (Taper Point)
SA 2030	2-0	30"	76 cm	30mm	1/2 Circle Round Body (Taper Point)
SA 2035	2-0	30"	76 cm	35mm	1/2 Circle Round Body (Taper Point)
SSA 2040	2-0	30"	76 cm	40mm	Straight Round Body
SSA 2050	2-0	30"	76 cm	50mm	Straight Round Body
SB 2013	2-0	30"	76 cm	13mm	3/8 Circle Round Body(Taper Point)
SB 2016	2-0	30"	76 cm	16mm	3/8 Circle Round Body(Taper Point)
SB 2020	2-0	30"	76 cm	20mm	3/8 Circle Round Body(Taper Point)
SB 2030	2-0	30"	76 cm	30mm	3/8 Circle Round Body(Taper Point)
SB 2045	2-0	30"	76 cm	45mm	3/8 Circle Round Body(Taper Point)
SC 2013	2-0	30"	76 cm	13mm	3/8 Circle Cutting
SC 2016	2-0	30"	76 cm	16mm	3/8 Circle Cutting
SC 2016DA	2-0	30"	76 cm	16mm	3/8 Circle Cutting Double Arm
SC 2022	2-0	30"	76 cm	22mm	3/8 Circle Cutting
SC 2026	2-0	30"	76 cm	26mm	3/8 Circle Cutting
SC 2045	2-0	30"	76 cm	45mm	3/8 Circle Cutting
SSC 2040	2-0	30"/36"	76/ 90cm	40mm	Straight Cutting
SSC 2050	2-0	30"/36"	76/ 90cm	50mm	Straight Cutting
SSC 2060	2-0	30"/36"	76/ 90cm	60mm	Straight Cutting
SD 2016	2-0	30"	76 cm	16mm	3/8 Circle Reverse Cutting
SD 2026	2-0	30"	76 cm	26mm	3/8 Circle Reverse Cutting
SD 2045	2-0	30"	76 cm	45mm	3/8 Circle Reverse Cutting
SD 2050	2-0	30"	76 cm	50mm	3/8 Circle Reverse Cutting
SE 2017	2-0	30"	76 cm	17mm	1/2 or 3/8 Circle Taper Cut

SMB Code	U.S.P. Size	Length Inches	Length Metric	Needle Length	Needle Curvature and Profile
SE 2025	2-0	30"	76 cm	25mm	1/2 or 3/8 Circle Taper Cut
SA 1020	0	30"	76cm	20mm	1/2 Circle Round Body (Taper Point)
SA 1025	0	30"	76cm	25mm	1/2 Circle Round Body (Taper Point)
SA 1026	0	30"	76cm	26mm	1/2 Circle Round Body (Taper Point)
SA 1030	0	30"	76cm	30mm	1/2 Circle Round Body (Taper Point)
SA 1035	0	30"	76cm	35mm	1/2 Circle Round Body (Taper Point)
SB 1016	0	30"	76cm	16mm	3/8 Circle Round Body (Taper Point)
SB 1020	0	30"	76cm	20mm	3/8 Circle Round Body (Taper Point)
SB 1030	0	30"	76cm	30mm	3/8 Circle Round Body (Taper Point)
SB 1045	0	30"	76cm	45mm	3/8 Circle Round Body (Taper Point)
SE 1025	0	30"	76cm	25mm	1/2 or 3/8 Circle Taper Cut
SSC 1040	0	30"/ 36"	76/90cm	40mm	Straight Cutting
SSC 1050	0	30"	76cm	50mm	Straight Cutting
SA 0120	1	30"	76cm	20mm	1/2 Circle Round Body (Taper Point)
SA 0125	1	30"	76cm	25mm	1/2 Circle Round Body (Taper Point)
SA 0126	1	30"	76cm	26mm	1/2 Circle Round Body (Taper Point)
SA 0130	1	30"	76cm	30mm	1/2 Circle Round Body (Taper Point)
SA 0135	1	30"	76cm	35mm	1/2 Circle Round Body (Taper Point)
SB 0120	1	30"	76cm	20mm	3/8 Circle Round Body (Taper Point)
SB 0130	1	30"	76cm	30mm	3/8 Circle Round Body (Taper Point)
SC 0160	1	30"	76cm	60mm	3/8 Circle Cutting
SD 0145	1	30"	76cm	45mm	3/8 Circle Reverse Cutting
SD 0150	1	30"	76cm	50mm	3/8 Circle Reverse Cutting
SE 0125	1	30"	76cm	25mm	1/2 or 3/8 Circle Taper Cut
SSC 0140	1	30"/36"	76/ 90cm	40mm	Straight Cutting
SSC 0150	1	30"/36"	76/ 90cm	50mm	Straight Cutting