



SMB Code	U.S.P. Size	Length Inches	Length Metric	Needle Length	Needle Curvature and Profile
PB 6010	6-0	30"/36"	76/90cm	10mm	3/8 Circle Round Body (Taper Point)
PB 6012	6-0	30"/36"	76/90cm	12mm	3/8 Circle Round Body (Taper Point)
PB 5010	5-0	30"/36"	76/90cm	10mm	3/8 Circle Round Body (Taper Point)
PB 5010 DA	5-0	30"/36"	76/90cm	10mm	3/8 Circle Round Body (Taper Point) Double Arm
PB 5012	5-0	30"/36"	76/90cm	12mm	3/8 Circle Round Body (Taper Point)
PB 5012 DA	5-0	30"/36"	76/90cm	12mm	3/8 Circle Round Body (Taper Point) Double Arm
PB 5016	5-0	30"/36"	76/90cm	16mm	3/8 Circle Round Body (Taper Point)
PC 5016	5-0	30"	76cm	16mm	3/8 Circle Cutting
PE 5017	5-0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut
PE 5017 DA	5-0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut Double Arm
PA 4016	4-0	30"/36"	76/90cm	16mm	1/2 Circle Round Body (Taper Point)
PA 4025	4-0	30"/36"	76/90cm	25mm	1/2 Circle Round Body (Taper Point)
PB 4010	4-0	30"/36"	76/90cm	10mm	3/8 Circle Round Body (Taper Point)
PB 4010 DA	4-0	30"/36"	76/90cm	10mm	3/8 Circle Round Body (Taper Point) Double Arm
PB 4012	4-0	30"/36"	76/90cm	12mm	3/8 Circle Round Body (Taper Point)
PB 4012 DA	4-0	30"/36"	76/90cm	12mm	3/8 Circle Round Body (Taper Point) Double Arm
PB 4016	4-0	30"/36"	76/90cm	16mm	3/8 Circle Round Body (Taper Point)
PC 4016	4-0	30"	76cm	16mm	3/8 Circle Cutting
PE 4017	4-0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut
PE 4017 DA	4-0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut Double Arm
PE 4045	4-0	30"/36"	76/90cm	45mm	1/2 or 3/8 Circle Taper Cut
PA 3016	3-0	30"	76cm	16mm	1/2 Circle Round Body (Taper Point)
PA 3016 DA	3-0	30"/36"	76/90cm	16mm	1/2 Circle Round Body (Taper Point) Double Arm
PA 3025	3-0	30"/36"	76/90cm	25mm	1/2 Circle Round Body (Taper Point)
PA 3025 DA	3-0	30"/36"	76/90cm	25mm	1/2 Circle Round Body (Taper Point) Double Arm
PA 3030	3-0	30"	76cm	30mm	1/2 Circle Round Body (Taper Point)
PB 3012	3-0	30"/36"	76/90cm	12mm	3/8 Circle Round Body (TaperPoint)
PB 3016	3-0	30"/36"	76/90cm	16mm	3/8 Circle Round Body (TaperPoint)
PC 3016	3-0	30"/36"	76/90cm	16mm	3/8 Circle Cutting
PE 3017	3-0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut
PE 3017 DA	3-0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut Double Arm
PE 3025	3-0	30"/36"	76/90cm	25mm	1/2 or 3/8 Circle Taper Cut
PE 3025 DA	3-0	30"/36"	76/90cm	25mm	1/2 or 3/8 Circle Taper Cut Double Arm
PA 2025	2-0	30"/36"	76/90cm	25mm	1/2 Circle Round Body (Taper Point)
PA 2030	2-0	30"	76cm	30mm	1/2 Circle Round Body (Taper Point)
PB 2016	2-0	30"/36"	76/90cm	16mm	3/8 Circle Round Body (TaperPoint)
PC 2016	2-0	30"	76cm	16mm	3/8 Circle Cutting
PE 2017	2-0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut
PE 2017 DA	2-0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut Double Arm



THE DEPENDABLE SUTURE  
TEEVRA CHAUKAS



Green/White Braided  
Silicon Coated Polyester

SMB Code	U.S.P. Size	Length Inches	Length Metric	Needle Length	Needle Curvature and Profile
PE 2025	2-0	30"/36"	76/90cm	25mm	1/2 or 3/8 Circle Taper Cut
PE 2025 DA	2-0	30"/36"	76/90cm	25mm	1/2 or 3/8 Circle Taper Cut Double Arm
PE 2045	2-0	30"/36"	76/90cm	45mm	1/2 or 3/8 Circle Taper Cut
PA 1025	0	30"/36"	76/90cm	25mm	1/2 Circle Round Body (Taper Point)
PA 1030	0	30"	76cm	30mm	1/2 Circle Round Body (Taper Point)
PA 1030DA	0	30"/36"	76/90cm	30mm	1/2 Circle Round Body (Taper Point) Double Arm
PE 1017	0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut
PE 1025	0	30"/36"	76/90cm	25mm	1/2 or 3/8 Circle Taper Cut
PE 1045	0	30"	76cm	45mm	1/2 or 3/8 Circle Taper Cut
PA 1025	0	30"/36"	76/90cm	25mm	1/2 Circle Round Body (Taper Point)
PA 0130DA	1	30"	76cm	30mm	1/2 Circle Round Body (Taper Point) Double Arm
PE 0125	1	30"/36"	76/90cm	25mm	1/2 or 3/8 Circle Taper Cut
PE 0145	1	30"	76cm	45mm	1/2 or 3/8 Circle Taper Cut
PE 0245	2	30"	76cm	45mm	1/2 or 3/8 Circle Taper Cut



THE DEPENDABLE SUTURE  
TEEVRA CHAUKAS



Green Braided PTFE  
Coated Polyester

SMB Code	U.S.P. Size	Length Inches	Length Metric	Needle Length	Needle Curvature and Profile
PTB 6010	6-0	30"/36"	76/90cm	10mm	3/8 Circle Round Body (Taper Point)
PTB 6010DA	6-0	30"/36"	76/90cm	10mm	3/8 Circle Round Body (Taper Point) Double Arm
PTB 5012	5-0	30"/36"	76/90cm	12mm	3/8 Circle Round Body (Taper Point)
PTB 5012DA	5-0	30"/36"	76/90cm	12mm	3/8 Circle Round Body (Taper Point) Double Arm
PTB 5013	5-0	30"/36"	76/90cm	13mm	3/8 Circle Round Body (Taper Point)
PTB 5013DA	5-0	30"/36"	76/90cm	13mm	3/8 Circle Round Body (Taper Point) Double Arm
PTB 5016	5-0	30"/36"	76/90cm	16mm	3/8 Circle Round Body (Taper Point)
PTB 5020	5-0	30"/36"	76/90cm	20mm	3/8 Circle Round Body (Taper Point)
PTB 5020DA	5-0	30"/36"	76/90cm	20mm	3/8 Circle Round Body (Taper Point) Double Arm
PTC 5016	5-0	30"	76cm	16mm	3/8 Circle Cutting
PTE 5017	5-0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut
PTE 5017DA	5-0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut Double Arm
PTA 4025	4-0	30"/36"	76/90cm	25mm	1/2 Circle Round Body (Taper Point)
PTB 4012	4-0	30"/36"	76/90cm	12mm	3/8 Circle Round Body (Taper Point)
PTB 4012DA	4-0	30"/36"	76/90cm	12mm	3/8 Circle Round Body (Taper Point) Double Arm
PTB 4013	4-0	30"/36"	76/90cm	13mm	3/8 Circle Round Body (Taper Point)
PTB 4013DA	4-0	30"/36"	76/90cm	13mm	3/8 Circle Round Body (Taper Point) Double Arm
PTB 4016	4-0	30"/36"	76/90cm	16mm	3/8 Circle Round Body (Taper Point)
PTB 4020	4-0	30"/36"	76/90cm	20mm	3/8 Circle Round Body (Taper Point)



SMB Code	U.S.P. Size	Length Inches	Length Metric	Needle Length	Needle Curvature and Profile
PTC 4016	4-0	30"	76cm	16mm	3/8 Circle Cutting
PTE 4017	4-0	30"	76cm	17mm	1/2 or 3/8 Circle Taper Cut
PTE 4017DA	4-0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut Double Arm
PTA 3025	3-0	30"/36"	76/90cm	25mm	1/2 Circle Round Body (Taper Point)
PTA 3025DA	3-0	30"/36"	76/90cm	25mm	1/2 Circle Round Body (Taper Point) Double Arm
PTB 3012	3-0	30"/36"	76/90cm	12mm	3/8 Circle Round Body (Taper Point)
PTB 3016	3-0	30"/36"	76/90cm	16mm	3/8 Circle Round Body (Taper Point)
PTB 3020	3-0	30"/36"	76/90cm	20mm	3/8 Circle Round Body (Taper Point)
PTC 3016	3-0	30"	76cm	16mm	3/8 Circle Cutting
PTE 3017DA	3-0	30"/36"	75/76/90cm	17mm	1/2 or 3/8 Circle Taper Cut Double Arm
PTE 3025	3-0	30"/36"	75/76/90cm	25mm	1/2 or 3/8 Circle Taper Cut
PTE 3025 DA	3-0	30"/36"	75/76/90cm	25mm	1/2 or 3/8 Circle Taper Cut Double Arm
PTA 2025	2-0	30"/36"	76/90cm	25mm	1/2 Circle Round Body (Taper Point)
PTA 2030	2-0	30"/36"	76/90cm	30mm	1/2 Circle Round Body (Taper Point)
PTB 2016	2-0	30"/36"	76/90cm	16mm	3/8 Circle Round Body (Taper Point)
PTB 2020	2-0	30"/36"	76/90cm	20mm	3/8 Circle Round Body (Taper Point)
PTC 2016	2-0	30"	76cm	16mm	3/8 Circle Cutting
PTE 2017	2-0	30"/36"	76/90cm	17mm	1/2 or 3/8 Circle Taper Cut
PTE 2017DA	2-0	30"/36"	75/76/90cm	17mm	1/2 or 3/8 Circle Taper Cut Double Arm
PTE 2025	2-0	30"/36"	75/76/90cm	25mm	1/2 or 3/8 Circle Taper Cut
PTE 2025DA	2-0	30"/36"	75/76/90cm	25mm	1/2 or 3/8 Circle Taper Cut Double Arm
PTA 1025	0	30"/36"	76/90cm	25mm	1/2 Circle Round Body (Taper Point)
PTA 1030	0	30"/36"	76/90cm	30mm	1/2 Circle Round Body (Taper Point)
PTA 1030DA	0	30"/36"	76/90cm	30mm	1/2 Circle Round Body (Taper Point) Double Arm
PTE 1045	0	30"/36"	76/90cm	45mm	1/2 or 3/8 Circle Taper Cut
PTA 0130	1	30"/36"	76/90cm	30mm	1/2 Circle Round Body (Taper Point)
PTE 0145	1	30"/36"	76/90cm	45mm	1/2 or 3/8 Circle Taper Cut
PTE 0245	2	30"/36"	76/90cm	45mm	1/2 or 3/8 Circle Taper Cut
PTE 0555	5	30"	75cm	55mm	1/2 or 3/8 Circle Taper Cut