



SYRINGE
HEALING LIFE

SMB[®] Luer Lock Disposable

**Sterile, Single use, Hypodermic Syringe
with or without Needle**

- Easy to draw up viscous allergenic extracts
- Only 0.0021 ml dead space
- 98% less antigen waste
- Smooth Plunger
- Clear Barrel
- Easy to read scale markings
- Ultra sharp needle



Benefits of Luer Lock Syringe

- The luer lock helps to ensure that the needle or tube attached will remain secure throughout the administration process of a treatment.
- Securing a proper connection between needle and patient, to ensure that there are no issues when it comes to the venous lines or ports to the arteriolar lines.
- Luer lock on your syringe, avoids the chances the needle could be lost in the patient.
- This technique avoids contamination occurring as the port can be cleaned and fresh syringes can be used instead of injecting patients with a port that has been used multiple times with different patients.
- The luer lock ensures a tight fit between the connectors and the needles.



Unit I
13, 33-36, Prem Industrial Estate,
Subhash Road, Jogeshwari (E),
Mumbai 400 060.

Unit II
Plot No.156 GIDC Umbergaon, 396170, Gujarat, INDIA.

Tel.: +91-22-2824 6922 / 2834 5637
E-Mail: sales@smbcorp.com | www.smbcorp.com