


"FOR A BETTER MOTHERHOOD, SAY YES TO IUDS, NO TO ABORTION"

SMB has its Global Presence In 140 countries
spanning all Continents with more than
65 million satisfied IUD users.



A DNV GL COMPANY



Indicative selection criteria of the Models Standard & Slek

Size	Vertical Arm Length	Horizontal Width	Approx. Sound Measuring Range
Standard	35.5 mm	19 mm	6-9 cm
Slek	29.5 mm	19 mm	5-8 cm

ADVANTAGES :

Proven, Safe, Reversible Contraception | Easy to load | Long Intrauterine life | Hormone/Steroid free | Highly convenient | Immediately effective, can be used as emergency contraceptive | Daily attention not required | Safe contraceptive for breast feeding women | No touch technique for insertion | Can be used by women who cannot use Steroid Contraceptives | Safe for all women of child bearing age | Safe for use after abortion or miscarriage.