

TCu 380A9


Long acting Contraceptive for **5 years**
NORMAL/MINI/MAXI

"FOR A BETTER MOTHERHOOD, SAY YES TO IUDS, NO TO ABORTION"

SMB has its Global Presence In **140** countries
 spanning all Continents with more than
65 million satisfied IUD users.



A DNV GL COMPANY



Indicative Selection Criteria of the Models Mini, Normal, Maxi

Size	Vertical Length	Horizontal Width	Approx. Sound Measuring Range
Mini	30.5 mm	36 mm	6 - 7.5 cm
Normal	33 mm	45 mm	7 - 8.5 cm
Maxi	38 mm	53 mm	8 - 9 cm

ADVANTAGES :

Proven, Safe, Reversible Contraception | Easy to load | Long Intrauterine life | Hormone/Steroid free | Highly convenient | Immediately effective, can be used as emergency contraceptive | Daily attention not required | Safe contraceptive for breast feeding women | No touch technique for insertion | Can be used by women who cannot use Steroid Contraceptives | Safe for all women of child bearing age | Safe for use after abortion or miscarriage.